


Product Datasheet

Human Fc Protein (orb334902)

Description	Crystallizable fragments composed of the carboxy-terminal halves of both IMMUNOGLOBULIN HEAVY CHAINS...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IgG1 Fc
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-6×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 26.9 kDa. The protein migrates as 30-33 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Glu 99 - Lys 330
Purity	95%
Source	Human IgG Fc, His Tag (orb334902) is expressed from human 293 cells (HEK293). It contains AA Glu 99 - Lys 330 (Accession # P01857-1).
MW	26.9 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human IgG
Fc protei...