


Product Datasheet

Human IL-17F (H161R) Protein (orb334906)

Description

Human IL-17F (H161R), His Tag (ILF-H4240) is



contains AA...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

IL-17F

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag

N-6×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 16.8 kDa. The protein migrates as 23 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

● -20°C to -70°C for 12 months in lyophilized state; ● -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Source


Human IL-17F (H161R), His Tag (orb334906) is expressed from human 293 cells (HEK293). It contains AA Arg 31 - Gln 163 (Accession # Q96PD4-1 (H161R)).

MW


16.8 kDa

Uniprot ID
Q96PD4
NCBI
AAH70124.1
Expiration Date

6 months from date of receipt.



SDS-PAGE analysis of Human IL-17F protei...



Loaded Human IL-17 RA, Fc Tag on Protein...