

Product Datasheet

Mouse IL-4 R alpha / CD124 Protein (orb334907)

Description Mouse IL-4 R alpha, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Ile 26 - Arg 2...

Reactivity Mouse

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target IL-4 R alpha

Preservatives Tris with Glycine, Arginine and NaCl, pH7.5

Form/Appearance Powder

Storage -20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag C-Fc

Note For research use only

Application notes This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 51 kDa. The protein migrates as 60-66 kDa under reducing (R) condition (SDS-PAGE).

Purity 95%


Source Mouse IL-4 R alpha, Fc Tag (orb334907) is expressed from human 293 cells (HEK293). It contains AA Ile 26 - Arg 233 (Accession # NP_001008700).

MW 51 kDa


Uniprot ID **P16382**

NCBI **NP_001008700**

Expiration Date 6 months from date of receipt.



SDS-PAGE analysis of Mouse IL-4 R alpha ...



Immobilized Mouse IL-4, His Tag (Cat. No...