


Product Datasheet

Mouse CD40 / TNFRSF5 Protein (orb334935)

Description	Mouse CD40, His Tag is expressed from human 293 cells (HEK293). It contains AA Val 24 - Arg 193 (Acc...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	CD40
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 20.8 kDa. The protein migrates as 25-30 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	P27512
Purity	95%
Source	Mouse CD40, His Tag (orb334935) is expressed from human 293 cells (HEK293). It contains AA Val 24 - Arg 193 (Accession # P27512-1).
MW	20.8 kDa
Uniprot ID	P27512
NCBI	P27512
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human
CD40
protein