


Product Datasheet

Mouse TIGIT Protein (orb334938)

Description	T-cell immunoreceptor with Ig and ITIM domains (TIGIT) is also known as V-set and immunoglobulin dom...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	TIGIT
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 39.5 kDa. The protein migrates as 45-55 kDa under reducing (R) condition (SDS-PAGE).
Protein Sequence	NP_001139797.1
Purity	95%
Source	Mouse TIGIT, Fc Tag (orb334938) is expressed from human 293 cells (HEK293). It contains AA Gly 26 - Thr 143 (Accession # NP_001139797.1).
MW	39.5 kDa
NCBI	NP_001139797.1
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human TIGIT
protein