


Product Datasheet

VEGFR2 antibody (orb334966)

Description	Goat polyclonal to VEGFR2. VEGFR2 known as kinase insert domain receptor, is a type III receptor tyr...
Species/Host	Goat
Reactivity	Bovine, Canine, Donkey, Equine, Feline, Gallus, Goat, Guinea pig, Hamster, Human, Mouse, Other, Porcine, Primate, Rabbit, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Purified recombinant peptide derived from the N-terminus (residues 20-125 aa) of human VEGFR2 produced in E. coli.
Target	VEGFR2
Preservatives	PBS, 20% glycerol and 0.05% sodium azide
Form/Appearance	Polyclonal antibody supplied as a 200 µl (3 mg/ml) aliquot in PBS, 20% glycerol and 0.05% sodium azide. This antibody is epitope-affinity purified from goat antiserum.
Concentration	3 mg/ml
Storage	Store at -20°C for long-term storage. Store at 2-8°C for up to one month.
Note	For research use only
Application notes	The antibody solution should be gently mixed before use.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2,000, IF:1:25-1:250
Expiration Date	12 months from date of receipt.



Western blot analysis using VEGFR2 antibody...