


Product Datasheet

INS antibody (orb334979)

Description	Goat polyclonal to Insulin. Insulin is a 6 kDa peptide, first synthesized as a precursor molecule, p...
Species/Host	Goat
Reactivity	Bovine, Canine, Donkey, Equine, Feline, Gallus, Goat, Guinea pig, Hamster, Human, Mouse, Other, Porcine, Primate, Rabbit, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	IHC-Fr, IHC-P, WB
Immunogen	Purified recombinant human Insulin produced in E. coli as a fusion protein.
Target	INS
Preservatives	PBS, 20% glycerol and 0.05% sodium azide
Form/Appearance	Polyclonal antibody supplied as a 200 µl (3 mg/ml) aliquot in PBS, 20% glycerol and 0.05% sodium azide. This antibody is epitope-affinity purified from goat antiserum.
Concentration	3 mg/ml
Storage	Store at -20°C for long-term storage. Store at 2-8°C for up to one month.
Note	For research use only
Application notes	The antibody solution should be gently mixed before use.
Protein Sequence	MFVNQHLCGSHLVEALYLVCGERGFFYTPKTGIVEQCCTSICSLYQLENYCN
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB:1:250-1:2,000, IHC-P:1:250-1:1,000, IHC-F:1:250-1:1,000
Expiration Date	12 months from date of receipt.



Western blot analysis of staining of MBP...