



Product Datasheet

RFP antibody (orb334992)


Description	Goat polyclonal antibody to RFP (Red fluorescent protein). RFP is derived from DsRed, an engineered ...
Species/Host	Goat
Conjugation	Unconjugated
Tested Applications	IEM, IF, IHC-Fr, IHC-P, WB
Immunogen	Purified recombinant RFP peptide produced in E. coli.
Target	RFP
Preservatives	PBS, 20% glycerol and 0.05% sodium azide
Form/Appearance	Polyclonal antibody supplied as a 100 µl (3 mg/ml) aliquot in PBS, 20% glycerol and 0.05% sodium azide. This antibody is epitope-affinity purified from goat antiserum.
Concentration	3 mg/ml
Storage	Store at -20°C for long-term storage. Store at 2-8°C for up to one month.
Note	For research use only
Application notes	The antibody solution should be gently mixed before use.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2,000, IF:1:50-1:500, IHC-P:1:50-1:500, IHC-F:1:50-1:500, IEM:1:50-1:500
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of staining of WB2...



Immunofluorescence analysis of 293HEK ce...