



Product Datasheet

Venus antibody (orb334993)


Description	Goat polyclonal antibody to Venus. Venus is a genetic mutant of green fluorescent protein (GFP), der...
Species/Host	Goat
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	IEM, IF, IHC-Fr, IHC-P, WB
Immunogen	Purified recombinant GFP peptide produced in E. coli.
Target	Venus
Preservatives	PBS, 20% glycerol and 0.05% sodium azide
Form/Appearance	Antibody supplied as a 100 µl aliquot in PBS, 20% glycerol and 0.05% sodium azide. This antibody is epitope-affinity purified from goat antiserum.
Concentration	3 mg/ml
Storage	Store at -20°C for long-term storage. Store at 2-8°C for up to one month.
Note	For research use only
Application notes	The antibody solution should be gently mixed before use.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2,000, IF:1:50-1:1,000, IHC-P:1:50-1:1,000, IHC-F:1:50-1:1,000, IEM:1:50-1:1,000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of staining of IF ...



Immunofluorescence analysis of 293HEK ce...