

## Product Datasheet

### Prostaglandin E2 receptor antibody (FITC) (orb3350)

**Description**

Rabbit polyclonal antibody against

FITC

**Species/Host**

Rabbit

**Reactivity**

Human, Mouse, Rat

**Conjugation**

FITC

**Tested**

IF

**Applications**
**Immunogen**

KLH conjugated synthetic peptide  
derived from human  
EP4R/Prostaglandin E Receptor EP4

**Target**

PTGER4

**Preservatives**

0.01M TBS(pH7.4) with 1% rAlbumin,  
0.03% Proclin300 and 50% Glycerol.

**Form/Appearance**

Lyophilized or Liquid

**Concentration**

- 100 µg (in 100 µl): 1 mg/ml

**Storage**

Store at -20 °C for one year. Avoid  
repeated freeze/thaw cycles. The  
lyophilized antibody is stable at room  
temperature for at least one month and  
for greater than a year when kept at -  
20°C. When reconstituted in sterile pH  
7.4 0.01M PBS or diluent of antibody  
the antibody is stable for at least two  
weeks at 2-4 °C.

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Polyclonal

**MW**

53kDa

**Uniprot ID**
**P35408**
**Expiration Date**

12 months from date of receipt.