


Product Datasheet

PLCE1 antibody (orb335375)

Description	Rabbit polyclonal antibody to PLCE1 also known as PLC epsilon1 and Phospholipase C epsilon1, whose e...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-2000, IHC-P:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293 cells
US...