


Product Datasheet

LYSMD3 antibody (orb33543)

Description	LYSMD3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human LysM and putative peptidoglycan-binding domain-containing protein 3 protein (1-127AA)
Target	LYSMD3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7Z3D4
Dilution Range	WB:1:1000-5000
Expiration Date	12 months from date of receipt.



Western blot
analysis of
Mouse heart
ti...