



Product Datasheet

CTIP (Phospho-Ser327) antibody (orb335672)


Description	Rabbit polyclonal antibody to CTIP (Phospho-Ser327)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:1000-3000, IHC-P: 1:50-1:200, ELISA(peptide): 1:20000-1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis of H2O2 treated 29...



Immunohistochemical staining of Rat inte...



Immunohistochemical staining of Mouse br...