


Product Datasheet

PKA alpha CAT (Phospho-Ser338) antibody (orb335685)

Description	Rabbit polyclonal antibody to PKA alpha CAT (Phospho-Ser338)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:1000-3000
Expiration Date	12 months from date of receipt.



Western blot analysis of Insulin treated...