



Product Datasheet

GCN2 (Phospho-Thr899) antibody (orb335822)

Description	Rabbit polyclonal antibody to GCN2 (Phospho-Thr899)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:1000-3000, IHC-P: 1:50-1:200, ELISA: 1:20000-1:40000
Expiration Date	12 months from date of receipt.



IHC-P
image of
Rat lung
tissue using
GCN...



IHC-P
image of
Mouse
brain tissue
using ...

