


Product Datasheet

p70 S6 Kinase (Phospho-Ser427) antibody (orb335845)

Description	Rabbit polyclonal antibody to p70 S6 Kinase (Phospho-Ser427)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:1000-3000
Expiration Date	12 months from date of receipt.



Western
blot
analysis of
HeLa
whole
cell...