


## Product Datasheet

### c-Ret (Phospho-Tyr981) antibody (orb335869)

|                            |                                                      |
|----------------------------|------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to c-Ret (Phospho-Tyr981) |
| <b>Species/Host</b>        | Rabbit                                               |
| <b>Reactivity</b>          | Human, Mouse, Rat                                    |
| <b>Conjugation</b>         | Unconjugated                                         |
| <b>Tested Applications</b> | WB                                                   |
| <b>Note</b>                | For research use only                                |
| <b>Clonality</b>           | Polyclonal                                           |
| <b>Dilution Range</b>      | WB: 1:1000-3000                                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.                      |



Western blot analysis of EGF treated A43...