


## Product Datasheet

### DAG1 (Phospho-Tyr892) antibody (orb335884)

|                            |                                                     |
|----------------------------|-----------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to DAG1 (Phospho-Tyr892) |
| <b>Species/Host</b>        | Rabbit                                              |
| <b>Reactivity</b>          | Human, Mouse                                        |
| <b>Conjugation</b>         | Unconjugated                                        |
| <b>Tested Applications</b> | WB                                                  |
| <b>Note</b>                | For research use only                               |
| <b>Clonality</b>           | Polyclonal                                          |
| <b>Dilution Range</b>      | WB: 1:1000-3000                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                     |



Western blot analysis of UV treated HeLa...