


Product Datasheet

Leptin Receptor (Phospho-Tyr1141) antibody (orb335892)

Description	Rabbit polyclonal antibody to Leptin Receptor (Phospho-Tyr1141)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:1000-3000
Expiration Date	12 months from date of receipt.



Western blot analysis of PMA treated Jur...