


Product Datasheet

Cyclin B1 (Phospho-Ser128) antibody (orb335933)

Description	Rabbit polyclonal antibody to Cyclin B1 (Phospho-Ser128)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:1000-3000
Expiration Date	12 months from date of receipt.



Western blot analysis of COS7 whole cell...