


Product Datasheet

SREBP-2 (Phospho-Ser455) antibody (orb335988)

Description	Rabbit polyclonal antibody to SREBP-2 (Phospho-Ser455)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:1000-3000
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 whole cel...