


Product Datasheet

LRP6 (Phospho-Thr1479) antibody (orb336011)

Description	Rabbit polyclonal antibody to LRP6 (Phospho-Thr1479)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:1000-3000
Expiration Date	12 months from date of receipt.



Western blot analysis of MDA-MB-435 whole cell lysate.