




Product Datasheet

CD73 antibody (orb33603)

Description	Rabbit polyclonal antibody to CD73
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	NT5E
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	63368
Uniprot ID	P21589
NCBI	NP_001191742.1 , NP_002517.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 cell lysate...



IHC-P of human placenta tissue (CD73 ant...