


## Product Datasheet

### LAT (Phospho-Tyr255) antibody (orb336104)

|                            |                                                    |
|----------------------------|----------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to LAT (Phospho-Tyr255) |
| <b>Species/Host</b>        | Rabbit                                             |
| <b>Reactivity</b>          | Human, Mouse, Rat                                  |
| <b>Conjugation</b>         | Unconjugated                                       |
| <b>Tested Applications</b> | WB                                                 |
| <b>Note</b>                | For research use only                              |
| <b>Clonality</b>           | Polyclonal                                         |
| <b>Dilution Range</b>      | WB: 1:1000-3000                                    |
| <b>Expiration Date</b>     | 12 months from date of receipt.                    |



Western blot analysis of H2O2 treated 29...