

## Product Datasheet

### SOX2 antibody (orb33646)

## Description

SOX2 antibody

## Species/Host

Rabbit

## Reactivity

Human, Mouse, Rat

## Conjugation

Unconjugated

## Tested

IHC-P, WB

## Applications

## Immunogen

KLH conjugated synthetic peptide derived from human SOX2. Please contact us for the exact immunogen sequence. The peptide is available as orb12925.

## Target

SOX2

## Form/Appearance

10 mM PBS, 0.02% sodium azide

## Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

## Storage

Store at 2-8°C for up to 6 months. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

## Note

For research use only

## Isotype

IgG

## Clonality

Polyclonal

## Purity

Polyclonal antibodies are purified by peptide affinity chromatography

## MW

34 kDa

## Uniprot ID

[P48431](#)

## NCBI

[003097](#), [003106](#), [31](#), [11](#)

## Entrez



[6657](#)

## Dilution Range


WB: 1:100-500, IHC-P: 1:100-500

## Expiration Date


12 months from date of receipt.



Immunohistochemical staining of human li...



Immunohistochemical staining of human st...



Immunohistochemical staining of paraffin...