




Product Datasheet

CD90 antibody (orb33652)

Description	Rabbit Polyclonal antibody againsts CD90
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	THY1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	17935
Uniprot ID	P04216
NCBI	NP_006279.2
Dilution Range	WB: 1:1000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



IHC-P human hepatocarcinoma tissue (THY1...



Western blot analysis of THY1 antibody i...