


Product Datasheet

Histone H2A (acetyl K9) antibody (orb336522)

Description	Rabbit polyclonal antibody to Histone H2A (acetyl K9)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:1000-3000, IHC-P: 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa whole cell...