


Product Datasheet

Histone H3 (mono+di+tri methyl K14) antibody (orb336524)

Description	Rabbit polyclonal antibody to Histone H3 (mono+di+tri methyl K14)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:1000-3000
Expiration Date	12 months from date of receipt.



Western blot analysis of HUVEC whole cel...