


Product Datasheet

RSK1 p90 (Phospho-T359+S363) antibody (orb336540)

Description	Rabbit polyclonal antibody to RSK1 p90 (Phospho-T359+S363)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:1000-3000, IHC-P: 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of COLO205 whole c...