


Product Datasheet

DRP1 (Phospho-Ser637) antibody (orb336568)

Description	Rabbit polyclonal antibody to DRP1 (Phospho-Ser637)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	A synthesized peptide derived from human DNML1(Accession O00429), corresponding to amino acid residues around phosphorylated Ser637
Target	Phospho-DRP1 (Ser637) Antibody detects endogenous levels of DRP1 only when phosphorylated at Ser637
Form/Appearance	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol
Storage	Store at -20 °C. Stable for 12 months from date of receipt
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	The antibody is from purified rabbit serum by affinity purification via sequential chromatography on phospho-peptide and non-phospho-peptide affinity columns
MW	82KD
Dilution Range	WB 1:1000-3000, IHC 1:50-1:200, IF/ICC 1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 whole cell ...