


Product Datasheet

Acetyl-alpha-Tubulin (Lys40) antibody (orb336569)

Description	Rabbit polyclonal antibody to Acetyl-alpha-Tubulin (Lys40)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:1000-3013
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa whole cell...