


Product Datasheet

Beclin 1 (Phospho-Ser93/96) antibody (orb336570)

Description	Rabbit polyclonal antibody to Beclin 1 (Phospho-Ser93/96)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:1000-3014
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 whole cell...