



Product Datasheet

p-PERK (Thr 982) antibody (orb336657)


Description	This p-PERK Ab is generated from rabbits immunized with a synthetic peptide. It is polyclonal and un...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:1000-3000, IHC-P: 1:50-1:200, IF: 1:100-500, ELISA: 1:20000-1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 whole cell ...



IF analysis of HepG2 cells lysates using...



Immunofluorescence analysis of A549 cells ...