


## Product Datasheet

### Doublecortin (Phospho-Ser334) antibody (orb336659)

|                            |                                            |
|----------------------------|--------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to Doublecortin |
| <b>Species/Host</b>        | Rabbit                                     |
| <b>Reactivity</b>          | Human                                      |
| <b>Conjugation</b>         | Unconjugated                               |
| <b>Tested Applications</b> | IHC-P, WB                                  |
| <b>Note</b>                | For research use only                      |
| <b>Clonality</b>           | Polyclonal                                 |
| <b>Dilution Range</b>      | WB: 1:1000-3000, IHC-P: 1:200              |
| <b>Expiration Date</b>     | 12 months from date of receipt.            |



Western blot analysis of K562 whole cell...