


Product Datasheet

INSRR antibody (orb336720)

Description	Rabbit polyclonal antibody to INSRR
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:1000-3000, IHC-P: 1:200, IF/ICC: 1:200-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 whole cell...