

Product Datasheet

FAM19A1 antibody (orb33675)

Description	FAM19A1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Protein FAM19A1 protein (20-133AA)
Target	FAM19A1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7Z5A9
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.