

## Product Datasheet

### FAM19A1 antibody (HRP) (orb33676)

|                            |                                                                  |
|----------------------------|------------------------------------------------------------------|
| <b>Description</b>         | FAM19A1 antibody (HRP)                                           |
| <b>Species/Host</b>        | Rabbit                                                           |
| <b>Reactivity</b>          | Human                                                            |
| <b>Conjugation</b>         | HRP                                                              |
| <b>Tested Applications</b> | ELISA                                                            |
| <b>Immunogen</b>           | Recombinant Human Protein FAM19A1 protein (20-133AA)             |
| <b>Target</b>              | FAM19A1                                                          |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol   |
| <b>Form/Appearance</b>     | Liquid                                                           |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C.<br>Avoid repeated freeze. |
| <b>Note</b>                | For research use only                                            |
| <b>Isotype</b>             | IgG                                                              |
| <b>Clonality</b>           | Polyclonal                                                       |
| <b>Uniprot ID</b>          | <a href="#">Q7Z5A9</a>                                           |
| <b>Dilution Range</b>      | ELISA                                                            |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                  |