


Product Datasheet

ROR2 antibody (orb336914)

Description	Rabbit polyclonal antibody to ROR2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:1000-3000, IF: 1:200-1:500, ELISA: 1:20000-1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa whole cell...