

Product Datasheet

SPATA19 antibody (Biotin) (orb33699)

Description	SPATA19 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Spermatogenesis-associated protein 19, mitochondrial protein (25-167AA)
Target	SPATA19
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7Z5L4
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.