


Product Datasheet

XPO4 antibody (orb337091)

Description	Rabbit polyclonal antibody to XPO4
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:1000-3000, IHC-P: 1:200, IF/ICC: 1:200-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HuvEc whole cel...