

Product Datasheet

UBE2G2 antibody (orb33711)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	UBE2G2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	18566
Uniprot ID	P60604
NCBI	NP_001189418.1, NP_872630.1, NP_003334.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

