




Product Datasheet

PGP9.5 antibody (orb33715)

Description	Rabbit polyclonal		
PGP9.5			
Species/Host	Rabbit		
Reactivity	Human, Mouse, Rat		
Conjugation	Unconjugated		
Tested Applications	FC, IF, IHC-P, WB		
Immunogen	Synthetic Peptide		
Target	UCHL1		
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.		
Note	For research use only		Immunohistochemical analysis of formalin... Immunofluorescence analysis of NCI-H460 ...
Application notes	Protein InformationName	Ubiquitin-protein hydrolase involved both in the processing of ubiquitin precursors and of ubiquitinated proteins. This enzyme is a thiol protease that recognizes and hydrolyzes a peptide bond at the C	Found in neuronal cell bodies and processes throughout the neocortex (at protein level). Expressed in neurons and cells of the diffuse Note=About neuroendocrine