



Product Datasheet

Inhibin beta -C antibody (orb337364)


Description	Rabbit polyclonal antibody to Inhibin beta -C
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:1000-3000, IHC-P: 1:200, ELISA: 1:20000-1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 whole cell ...



Immunohistochemical staining of Human Me...



Western blot analysis of extracts from M...