


Product Datasheet

Sarcoglycan-beta antibody (orb337401)

Description	Rabbit polyclonal antibody to Sarcoglycan-beta
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:1000-3000, IHC-P: 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat whole ce...