


## Product Datasheet

### ADCK4 antibody (orb337430)

|                            |                                     |
|----------------------------|-------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to ADCK4 |
| <b>Species/Host</b>        | Rabbit                              |
| <b>Reactivity</b>          | Human, Mouse, Rat                   |
| <b>Conjugation</b>         | Unconjugated                        |
| <b>Tested Applications</b> | WB                                  |
| <b>Note</b>                | For research use only               |
| <b>Clonality</b>           | Polyclonal                          |
| <b>Dilution Range</b>      | WB: 1:1000-3000                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.     |



Western blot analysis of HuvEc whole cel...