


Product Datasheet

OSR1 (Phospho-Thr185) antibody (orb337538)

Description	Rabbit polyclonal antibody to OSR1 (Phospho-Thr185)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:1000-3000, IHC-P: 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 whole cell...