


## Product Datasheet

### PIF1 antibody (orb338052)

|                            |                                    |
|----------------------------|------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to PIF1 |
| <b>Species/Host</b>        | Rabbit                             |
| <b>Reactivity</b>          | Human, Mouse, Rat                  |
| <b>Conjugation</b>         | Unconjugated                       |
| <b>Tested Applications</b> | IHC-P, WB                          |
| <b>Note</b>                | For research use only              |
| <b>Clonality</b>           | Polyclonal                         |
| <b>Dilution Range</b>      | WB: 1:1000-3000, IHC-P: 1:200      |
| <b>Expiration Date</b>     | 12 months from date of receipt.    |



Western blot analysis of HepG2 whole cel...