


## Product Datasheet

### B4GALT7 antibody (orb338072)

|                            |                                       |
|----------------------------|---------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to B4GALT7 |
| <b>Species/Host</b>        | Rabbit                                |
| <b>Reactivity</b>          | Human, Mouse, Rat                     |
| <b>Conjugation</b>         | Unconjugated                          |
| <b>Tested Applications</b> | IHC-P, WB                             |
| <b>Note</b>                | For research use only                 |
| <b>Clonality</b>           | Polyclonal                            |
| <b>Dilution Range</b>      | WB: 1:1000-3000, IHC-P: 1:200         |
| <b>Expiration Date</b>     | 12 months from date of receipt.       |



Western blot analysis of A549 whole cell...