

Product Datasheet

C5orf24 antibody (Biotin) (orb33865)

Description	C5orf24 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human UPF0461 protein C5orf24 protein (7-188AA)
Target	C5orf24
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7Z6I8
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.