



Product Datasheet

FOXO1 (phospho-S256) antibody (orb338819)


Description	Rabbit polyclonal antibody to FOXO1
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IF, IH, IP, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S256 of human FOXO1 protein. The exact sequence is proprietary.
Target	FOXO1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q12778, G3V7R4, Q9R1E0
Entrez	84482, 2308, 56458
Dilution Range	WB: 1:500-1000, IHC-P: 1:100-200, IF/ICC: 1:100-500, IP: 1:10-100
Expiration Date	12 months from date of receipt.



Western blot analysis of FOXO1 (Phospho-S256) antibody.



Immunohistochemical analysis of FOXO1 (Phospho-S256) antibody.



Immunofluorescent analysis of FOXO1 (Phospho-S256) antibody.